

# Spherotech, Inc.


1840 Industrial Dr., Suite 270, Libertyville, IL 60048

## SPHERO™ Ultra Rainbow Calibration Particles


A new set of calibration particles, **Ultra Rainbow Calibration Particles**, is now available for performance tracking of flow cytometers with fluorescent channels in the Far Red. The **Ultra Rainbow Calibration Particles** are available in either 6 peaks, 3.8µm or 5.05µm. The **Ultra Rainbow Calibration Particles** have improved resolution in the PE-Cy7, APC, and APC-Cy7 channels. These particles are supplied as a kit. One 3mL bottle of blank particles and a 2 mL bottle of fluorescent particles.

Product	Size, µm	Catalog No.	Unit
Ultra Rainbow Calibration, 6 peaks, 10 <sup>7</sup> /mL	3.5-3.9	URCP-38-2K	2 mL
Ultra Rainbow Calibration, 6 peaks, 10 <sup>7</sup> /mL	5.0-5.9	URCP-50-2K	2 mL

### Histograms of URCP-38-2K



### Histograms of URCP-50-2K



Tel.: 800-368 0822 or 847-680 8922; Fax: 847-680 8927; E-Mail: [service@spherotech.com](mailto:service@spherotech.com)

Visit us on the web at <http://www.spherotech.com>