




Spherotech

Let the Possibilities Flow

SPHERO™ Drop Delay Calibration Particles

Fine Tune Your Sorter with the Appropriate Attachment and Spherotech Drop Delay Calibration Particles


- Aids in the determination of the drop delay value
- Increases the accuracy and productivity while sorting
- Consists of 2mL at 1×10^8 particles/mL of a single population of fluorescent particles
- Contains a mixture of fluorophores which allow detection in multiple channel flow cytometers



Histogram of Drop Delay Calibration Particles
Catalog Number: DDCP-70-2



Initial Drop Delay Profile with Spherotech Drop Delay Calibration Particles



Optimized Drop Delay Profile with Spherotech Drop Delay Calibration Particles

To learn more about Spherotech Drop Delay Calibration Particles visit us online

www.spherotech.com

27845 Irma Lee Circle, Unit 101
Lake Forest, IL 60045-5100
(847) 680-8922 fax (847) 680-8927

